



**SPECIALIZED
MANUFACTURER
COMPANY**

IN MEDICAL EQUIPMENT INDUSTRY

**MEDITECH CORPORATION
HAS BEEN **ISO & CE CERTIFIED**
MEDICAL EQUIPMENT
MANUFACTURER COMPANY,
STRICTLY FOLLOWING **QMS,**
CRM AND CALL CENTRE
FACILITIES TO ENHANCE THE
CUSTOMER SATISFACTION.**

ISO 13485:2016 CERTIFIED PICL/QMS/0420/7466

ISO 23500-5:2019 PICL/QMFHRT/0420/4897

CE CERTIFICATE NO: PICL/CE/0420/3618

**ONE CALL SOLUTION
FACILITATES CONFIDENTIAL
AND TRANSPARENT SUPPORT.**

**QUALITY
MANAGEMENT**

CUSTOMER SATISFACTION

ISO STANDARDS



**MEDITECH CORPORATION
WORKS WITH PRIVATE,
PUBLIC & NOT-FOR-PROFIT
ORGANIZATIONS.**

CONTACT

MediTech Coporations
72/26-B2, Nadeppally
-Complex, Kochi -682026

Call Centre: 9846 746 764

www.meditechcorp.com
enquiry@meditechcorp.com



ECG MACHINE MTC E7100



9.7 INCH TILTED TOUCH SCREEN LCD DISPLAY
HIGH PERFORMANCE CHIP MODULE MAKES
ACCURATE ECG SIGNAL COLLECTION
12 LEAD SYNCHRONOUS ACQUISITION AND DISPLAY
SAMPLING RATE 32000Hz
INSTITUTION & PATIENT INFO ENTRY
BUILT IN 1000 ECG STORAGE, EXPANDABLE MEMORY
USB & SD CARD
REALTIME ECG FREEZING & REVIEW
FORMAT:-PDF/ PNG/ JPG/ GDT/ DICOM/XML/ HL7/
ZQ ECG
NETWORK :-DICOM, WIFI,LAN

ECG MACHINE MTC E617



COMPACT & EASY OPERATION
ACCURATE PACE PULSE IDENTIFICATION &
AUTOMATIC ANALYSIS FUNCTION
OPERATION:- 10 INCH GRAPHIC COLOR DISPLAY
(800*480), HIGH BATTERY BACKUP
MODES: - AUTO/MANUAL/R-R ANALYSIS/STORAGE
MODES
PRINT: -3Ch+++,6Ch,6Ch+,12Ch FORMAT
STORAGE : INBUILT 250 ECG STORAGE (REVIEW
& PRINT)
NETWORK :-LAN

ECG MACHINE MTC E317



5 INCH HIGH LCD DISPLAY WITH 8 HOURS BATTERY
BACKUP
ACCURATE ECG SIGNAL COLLECTION
SAMPLING RATE 32000Hz
INSTITUTION & PATIENT INFO ENTRY
BUILT IN 1000 ECG STORAGE, EXPANDABLE MEMORY
USB & SD CARD
FORMAT:-PDF/ PNG/ JPG/ GDT/ DICOM/XML/ HL7/
ZQ ECG
NETWORK :-DICOM, WIFI,LAN
EXTERNAL CONNECTIVITY:-PRINTER,KEYBOARD,MOUSE

PATIENT MONITOR MTC M910



ECG,RESP,SPO2,PR,NIBP,DUAL TEMP
12.1 INCH TFT COLORFUL AND CLEAR DISPLAY
4 HOURS BATTERY BACKUP
PATIENT DATA RECORD, DISPLAY MODES
THREE LEVEL AUDIBLE AND VISUAL ALARM
13 KIND OF ARRHYTHMIA ANALYSIS
15 KINDS OF DRUG DOSE CALCULATIONS
SUPPORT MASIMO SPO2-AUTHORIZATION PARTNER
NETWORK CONNECTIVITY

PATIENT MONITOR MTC M710



ECG, RESP, SPO2, PR, NIBP, DUAL TEMP
12.1 INCH WITH TRUE COLOR LCD DISPLAY
LARGE CAPACITY 4000mAh BATTERY
PLAYBACK AND BROWSE FUNTION FOR WAVEFORM
AND MONITOR
PATIENT DATA RECORD, DISPLAY MODES
AUTO DOUBLE ALARM WITH AUDIBLE AND VISIBLE
SIGNALS
NETWORK CONNECTIVITY

PATIENT MONITOR MTC M319



CONTINUOUSLY MEASURE SPO2, NIBP & PR
COMPACT AND PORTABLE
4.3" LCD DISPLAY, 6HOURS BACKUP
TRUE COLOR WIDE VIEWING CAPACITY
PATIENT DATA RECORD, DISPLAY MODES
BAR GRAPH & PLETHYSMOGRAPH DISPLAY
36 HOURS GRAPHICAL AND TABULAR TREND DATA
STORAGE

FETAL MONITOR MTC F372



TWINS MONITORING SYSTEM
PATIENT DATA ENTRY OPTION
10.4 INCH TFT LED DISPLAY WITH 2HOURS BATTERY
BACKUP
PATIENT DATA RECORD, DISPLAY MODES
DISPLAY & PRINT: FHR1,FHR2,TOCO,AFM
WAVEFORMS
CROSS CHANNEL VERIFICATION (CCV)
PDF EXPORT TO USB
12-CRYSTAL PULSE DOPPLER TRANSDUCER

INFUSION PUMP MTC I421

CONVINIENT, USER FRIENDLY, EASY TO OPERATE
USER FRIENDLY INSTALLATION PROCESS FOR I.V. SET
TUBE
INFUSED VOLUME DISPLAY
EASY PROGRAMMING, 4 HOURS BATTERY BACKUP
WORKS ON VARIOUS I.V.SETS
AIR DETECTOR AND PRESSURE SENSOR
PROCESS FOLLOWED BY AN IMMEDIATE ALARM
SIGNAL



SYRINGE PUMP MTC S213



SAFE, CONVINIENT, USER FRIENDLY, EASY TO OPERATE
LED DISPLAY FOR FLOW RATE, 4 HOURS BATTERY BACKUP
ACCEPT SYRINGES OF 10ml, 20ml, 30ml,AND 50/60ml
IDENFIES AND DISPLAY THE SYRINGE TYPE AUTOMATICALLY
MODES OF OPERATIONS- VOLUME/TIME,RATE,BODY WEIGHT
ACCURATE AND SENSITIVE PRESSURE SENSOR SHORTENS

DIALYSIS MACHINE



NIPRO IS ONE OF THE LARGEST MANUFACTURERS OF DIALYZERS IN THE WORLD, BUT UNDERSTANDS THAT THE DIALYSIS MACHINE IS THE HEART OF THE TREATMENT. THE SURDIAL™ 55PLUS IS AN EASY TO USE DIALYSIS MACHINE OFFERING PATIENTS THE BEST RENAL REPLACEMENT THERAPY TO INCREASE THEIR WELL BEING AND QUALITY OF LIFE.

DIALYSIS CHAIR

HIGHEST COMFORT, ATTRACTIVE DESIGN IN ADDITION TO FUNCTIONALITY AND SAFETY.

- INTELLIGENT CONTROL SYSTEM
- FULLY REMOTE CONTROL
- WATER & DUST PROOF REMOTE CONTROL



PORTABLE RO PLANTS



DIA-PURE SERIES PORTABLE RO SYSTEM WITH ULTRA FILTRATION PRE FILTERS IS A COMPLETE REVERSE OSMOSIS WATER TREATMENT SYSTEM DESIGNED TO OPERATE FOR TWO MACHINES WITH A TREATED WATER TANK. THE COMPLETE SYSTEM IS HOUSED WITHIN A COMPACT SOUND PROOFED FRAME PROVIDING QUIET OPERATIONS AND PORTABILITY.

BICARBONATE MIXER



DIA MIX SERIES BI CARBONATE MIXER DESIGNED TO PREPARE A LIQUID SODIUM BICARBONATE FORMULATION FOR THE PURPOSE OF MAKING DIALYSATE FOR HAEMODIALYSIS PATIENT THERAPIES. WE USED SS 304 AGITATOR.

WE PROVIDE CUSTOMIZATION OPTIONS FOR PRODUCTS TO SUIT THE NEEDS OF THE CUSTOMERS. IN ADDITION, WE OFFER STANDARD DESIGNS THAT ARE UNIVERSALLY USED AND TESTED IN SERVICE.

WASHING SCRUB



DIALYZER REPROCESSOR



DIALYZER REPROCESSING MACHINE ARE THE FIRST AUTOMATIC DIALYZER REPROCESSING MACHINE FIVE WORKSTATION. OUR PERFECTNESS COMES FROM PROFESSIONAL AND ADVANCED TECHNOLOGY, THAT MAKES OUR PRODUCTS LEGAL, SAFE, AND STABLE.

FULLY AUTOMATED RO PLANT



A COMPREHENSIVE FULLY AUTOMATED DIA SAFE SERIES WATER TREATMENT PLANT EXCLUSIVELY DESIGNED FOR DIALYSIS TREATMENT AS PER AAMI STANDARDS. THE SYSTEM MINIMIZES THE END-USER DAILY ACTIVITY AND NEGLIGENCE OF OPERATOR ISSUES. THE SYSTEM WILL MAINTAIN AUTOMATIC BACKWASH, RINSE AND REGENERATION OF SOFTNER UNITS, HENCE THE PERFORMANCE OF RO PLANT WILL BE MORE EFFICIENT AND ACCURATE, WHICH WILL INCREASE MEMBRANE LIFE. UF FILTER ULTRA FILTRATION IS USED FOR PRE-FILTRATION IN REVERSE –OSMOSIS PROCESS. ULTRA FILTRATION IS AN EFFECTIVE MEANS OF REDUCING PARTICULATES THAT CAN FOUL REVERSE OSMOSIS MEMBRANE.



ULTRA FILTER WATER TREATMENT PLANT (WATER TREATMENT PLANT)


WE ARE SPECIALIZED IN ULTRA FILTERED WATER TREATMENT PLANT MANUFACTURER IN KERALA, WE HAD INSTALLED LARGEST ULTRA FILTERED WATER TREATMENT SYSTEM, THE SYSTEM CONTAIN BASIC PRE FILTER LIKE SAND , IRON REMOVER, AND ACTIVATED CARBON FILTER THEN THE FILTERED WATER PASSES THROUGH ULTRA FILTERED MEMBRANE SUSPENDED SOLIDS.



SEWAGE/ EFFLUENT WATER TREATMENT PLANT (STP/ETP)

SEWAGE TREATMENT IS A TYPE OF WASTEWATER TREATMENT WHICH AIMS TO REMOVE CONTAMINANTS FROM SEWAGE TO PRODUCE AN EFFLUENT THAT IS SUITABLE FOR DISCHARGE TO THE SURROUNDING ENVIRONMENT OR AN INTENDED REUSE APPLICATION, THEREBY PREVENTING WATER POLLUTION FROM RAW SEWAGE DISCHARGES.

OUR INSTALLATIONS / CLIENTS IN KERALA



WHY MEDITECH?

- ✓ Competitive Price
- ✓ Turnkey Project Experience
- ✓ Qualified Engineers
- ✓ Largest installations in Kerala
- ✓ Largest installations in Private Hospitals
- ✓ Largest installations in Government Hospitals

MAJOR CLIENTS



CONTACT

MediTech Corporations
72/26-B2, Nadeppally
-Complex, Kochi -682026

www.meditechcorp.com

enquiry@meditechcorp.com

Call Centre: 9846 746 764

