



**SPECIALIZED
MANUFACTURER
COMPANY**

IN MEDICAL EQUIPMENT INDUSTRY

**MEDITECH CORPORATION
HAS BEEN **ISO & CE CERTIFIED**
MEDICAL EQUIPMENT
MANUFACTURER COMPANY,
STRICTLY FOLLOWING **QMS,**
CRM AND CALL CENTRE
FACILITIES TO ENHANCE THE
CUSTOMER SATISFACTION.**

ISO 13485:2016 CERTIFIED PICL/QMS/0420/7466

ISO 23500-5:2019 PICL/QMFHRT/0420/4897

CE CERTIFICATE NO: PICL/CE/0420/3618

**ONE CALL SOLUTION
FACILITATES CONFIDENTIAL
AND TRANSPARENT SUPPORT.**

**QUALITY
MANAGEMENT**

CUSTOMER SATISFACTION

ISO STANDARDS



**MEDITECH CORPORATION
WORKS WITH PRIVATE,
PUBLIC & NOT-FOR-PROFIT
ORGANIZATIONS.**

CONTACT

MediTech Coporations
72/26-B2, Nadeppally
-Complex, Kochi -682026

Call Centre: 9846 746 764

www.meditechcorp.com
enquiry@meditechcorp.com



ECG MACHINE

MTC E317



ABOUT:-

PATIENT DATA & HOSPITAL INFO ENTRY
SAMPLING RATE 3200Hz
PACE DETECTION
AC POWER & BATTERY
CONTINUOUS OPERATION TIME IN BATTERY
POWER : MORE THAN 8 HOURS

ECG REPORTS:-

AUTO ECG REPORT WITH & WITHOUT INTERPRETATION,
AUTO ECG RHYTHM LEAD ANALYSIS REPORT WITH &
WITHOUT INTERPRETATION, MANUAL MODE-REAL
TIME RHYTHM.

PRINT:-

BUILT IN PRINTING FORMAT: 3x4, 3x4+1R, 1x12.
EXTERNAL PRINTING FORMAT : 3x4, 3x4+1R, 6x2,
6x2+1R, 1x12.

SAVING OPTION:-

ECG SAVING OPTION: WITH OR WITHOUT PRINT
INBUILT 1000 ECG STORAGE. EXPANDABLE MEMORY
USB OR SD CARD. FORMAT: PDF/ PNG/ JPG/ GDT/
DICOM/XML/ HL7/ZQ ECG

DISPLAY:-

5 INCH HIGH RESOLUTION LCD DISPLAY

NETWORK:-

DICOM, LAN, WIFI

EXTRAS:-

HIGH PERFORMANCE CHIP MODULE MAKES ACCURATE-
ECG SIGNAL COLLECTION
KEYBOARD, MOUSE, EXTERNAL PRINTER STAND BY
MODE, AUTO SHUT DOWN

TECHNICAL SPECIFICATIONS

SAFETY INDEX:-

CLASS I, TYPE CF, WITH PROTECTIVE CIRCUIT OF DEFIBRILLATION
AND PACE-MAKING

POWER SPECIFICATIONS:-

AC POWER 100V-240V, 50Hz/60Hz $\pm$ 1Hz, 75VA
BATTERY RECHARGEABLE Li-ion BATTERY (10.8V/2200mAh)

APPEARANCE PARAMETERS:-

SIZE (L*W*H) 264.4mm*183.9mm*87.6mm
WEIGHT ABOUT 1.3 Kg

ENVIRONMENTAL CONDITIONS

OPERATION:-

AMBIENT TEMPERATURE	+5C -- +40C
AMBIENT HUMIDITY	20% -- 85% (NO CONDENSATION)
ATMOSPHERIC PRESSURE	57kpa -- 106kpa

SHIPMENT AND STORAGE:-

AMBIENT TEMPERATURE	-20 C-- +55C
AMBIENT HUMIDITY	10% -- 95%
ATMOSPHERIC PRESSURE	500hPa -- 1060hPa



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2, Nadeppally Complex
Kochi-682026
www.meditechcorp.com



ECG MACHINE

MTC E617



FEATURES:-

COMPACT & EASY OPERATION
ACCURATE PACE PULSE IDENTIFICATION &
AUTOMATIC ANALYSIS FUNCTION

OPERATION:-

AC POWER & BATTERY
HIGH CAPACITY 4400mAh BATTERY

DISPLAY:-

10 INCH GRAPHIC COLOR DISPLAY (800*480)
12 LEAD SIMULTANEOUS WAVEFORM DISPLAY

ECG REPORTS:-

AUTO/MANUAL/R-R ANALYSIS/STORAGE MODES
AUTOMATIC ANALYSIS FUNCTION
3Ch+++, 6Ch, 6Ch+, 12Ch FORMAT

SAVING OPTION:-

INBUILT 250 ECG STORAGE (REVIEW & PRINT)

NETWORK:-

LAN

TECHNICAL SPECIFICATIONS

SAFETY INDEX:-

CLASS I, TYPE CF, WITH ANTI-DEFIBRILLATION FUNCTION
ELECTROMAGNETIC COMPATIBILITY: CISPR 11, GROUP I, CLASS A

POWER SPECIFICATIONS:-

AC POWER 100V-240V, 50Hz/60Hz ± 1 Hz, 60VA
BATTERY RECHARGEABLE Li-ion BATTERY
(14.4V/2200mAh or 4400mAh)

APPEARANCE PARAMETERS:-

SIZE (L*W*H) 347 mm*293mm*83mm
WEIGHT ABOUT 4.8 Kg

ENVIRONMENTAL CONDITIONS

OPERATION:-

AMBIENT TEMPERATURE	+5C -- +40C
RELATIVE HUMIDITY	25% -- 80% (NO CONDENSATION)
ATMOSPHERIC PRESSURE	700hPa – 1060hPa

SHIPMENT AND STORAGE:-

AMBIENT TEMPERATURE	-20 C-- +55C
AMBIENT HUMIDITY	10% -- 95%
ATMOSPHERIC PRESSURE	500hPa – 1060hPa



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2, Nadeppally Complex
Kochi-682026
www.meditechcorp.com



FETAL MONITOR

MTC F372



MAIN FEATURES:-

MAIN BOARD WITH LINUX OS
12-CRYSTAL PULSE DOPPLER TRANSDUCER
EVENT MARKER:-ANNOTATE PERCEIVED FETAL-
MOVEMENTS OR ANY OTHER EVENTS
AUTOMATIC/MANUAL FETAL-
MOVEMENT DETECTION
PATIENT DATA ENTRY OPTION
HEART RATE OFFSET MODE : VISUALLY-
DISTINGUISH BETWEEN TWINS HEART RATE
CROSS CHANNEL VERIFICATION (CCV) FETAL OR-
MATERNAL HEART RATE
NST TIMER

OPERATION:-

AC POWER & BATTERY
LARGE CAPACITY 2200mAh Lithium BATTERY
(4-5 HOURS,WORKING TIME :2 HOURS)

DISPLAY:-

FOLDED 10.4 INCH HIGH BRIGHTNESS TFT LED
DISPLAY

MONITOR:-

FHR1,FHR2,TOCO,AFM WITH WAVEFORM
MODES: STANDARD INTERFACE AND BIG FONT
INTERFACE
FHR SIGNAL QUALITY INDICATOR HELPS OPTIMIZE
THE PROBE POSITION

NETWORK:-

LAN CONNECTIVITY

TECHNICAL SPECIFICATIONS

PRINT:-

FHR1,FHR2,TOCO,AFM
OFFSET OPTION:FHR2 FROM -30 TO +30
PAPER SPEED : 1cm/min, 2cm/min, 3cm/min
CAN SET WAVEFORM WIDTH SIZE:NORMAL,DEMI-BOLD,BOLD
PRINTER : THERMAL 112mm/ Z-fold (REAL TME & REVIEW)

SAVING OPTION:-

8640 HOURS MEMORY FOR FETAL TRACES
PDF EXPORT TO USB

EXTRAS:-

KEYBOARD,MOUSE

SAFETY INDEX:-

EQUIPMENT EMC : IEC60601- CLASS A OF THE GROUP I
CLASS II EQUIPMENT WITH THE ANTI-ELECTRIC SHOCK DEGREE
OF BF

POWER SPECIFICATIONS:-

AC POWER	100V-240V AC,50Hz/60Hz,75VA
ADAPTER	9V DC 5A
BATTERY	RECHARGEABLE Li-ion BATTERY/ 7.4V 2200mAh

APPEARANCE PARAMETERS:-

SIZE	36cm (W) *26cm (D) *8cm (H)
WEIGHT	ABOUT 2.8kg



Contact us now
+91 98467 46764

Meditech Corporations
72/26-B2,Nadeppally Complex
Kochi-682026
www.meditechcorp.com



INFUSION PUMP MTC I421



FEATURES:-

CONVINIENT, USER FRIENDLY, EASY TO OPERATE
USER FRIENDLY INSTALLATION PROCESS FOR-
I.V. SET TUBE
INFUSED VOLUME DISPLAY
VOICE NOTIFICATION FOR ALARM CONDITION
ERGONOMOMIC TOUCH PANEL KEYPAD FOR-
EASY OPERATIONS
COLOR POWER INDICATORS
WORKS ON VARIOUS I.V.SETS
AIR DETECTOR AND PRESSURE SENSOR-
PRECISELY DETECTS AIR-IN-LINE AND STOPS THE-
INFUSION
PROCESS FOLLOWED BY AN IMMEDIATE ALARM-
SIGNAL

OPERATION:-

AC POWER & BATTERY
BATTERY BACKUP : 4 HOURS

DISPLAY:-

LARGER LCD SCREEN FOR DISPLAY FACILITATES-
EASY PROGRAMMING
LCD DISPLAY WITH BACKLIGHT FOR MENU-
SETTINGS AND SYSTEM STATUS

TECHNICAL SPECIFICATIONS

APPEARANCE PARAMETERS:-

SIZE (L*H*D) 140mm*157mm*220mm
WEIGHT 1.8 Kg

POWER :-

AC POWER 100V -240 V,50/60 Hz, 30VA
BATTERY Li-ion BATTERY,1.6Ah,4 HOURS BACKUP
RATE 1ml/h -1,200ml/h
FLOW RATE ACCURACY WITHIN $\pm 5\%$
(AFTER CURRENT CALIBRATION)

OCCCLUSION:-

HIGH	800mmHg ± 200 mmHg (106.7kPa ± 26.7 kPa)
MEDIUM	500mmHg ± 100 mmHg (66.7kPa ± 13.3 kPa)
LOW	300mmHg ± 100 mmHg (40.7kPa ± 13.3 kPa)

ALARMS:-

INFUSION COMPLETE, OCCCLUSION, DOOR OPEN,AIR LINE-
WRONG SETTING,LOW BATTERY,NO AC POWER ETC.

AIR IN LINE DETECTOR :-

ULTRASONIC DETECTOR,SENSITIVITY > 25 micro l
MECHANICAL PRECISION WITHIN $\pm 2\%$
INFUSION VOLUME 1ml – 9999ml
KVO RATE 1ml/h – 5ml/h (1ml/h step)



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2,Nadeppally Complex
Kochi-682026
www.meditechcorp.com



PATIENT MONITOR

MTC M710



ABOUT:-

MEASURED PARAMETER: ECG, RESP, SPO2, PR, NIBP, DUAL TEMP
SIMPLE AND FRIENDLY OPERATING DISPLAY-INTERFACE
PLAYBACK AND BROWSE FUNCTION FOR WAVEFORM AND MONITOR PATIENT DATA
A FULL-SYNCHRONISTIC LEAD MULTI-CHANNEL ECG DISPLAY
REAL TIME ST SEGMENT ANALYSIS AND PACE-MAKER DETECTION MODE
DYNAMIC WAVEFORMS AND MULTI-LEAD ECG-WAVE-FORMS DISPLAY
PATIENT INFO INPUT MANAGEMENT FUNCTION
USB AND NETWORK PORT

OPERATION:-

AC POWER & BATTERY
LARGE CAPACITY 4000mAh BATTERY (MORE THAN 60 MINUTES BACKUP)

DISPLAY:-

12.1 INCH WITH TRUE COLOR LCD DISPLAY
WIDE VIEWING ANGLE
HIGH BRIGHTNESS LCD DISPLAY

MONITOR:-

DISPLAY MODES: STANDARD, TREND SCREEN, BIG FONT, OXY CRG
SUITABLE FOR ADULT, PEDIATRIC, NEONATE (BEDSIDE/OR/ICU/NICU/PICU)

DATA STORAGE:-

400 NIBP MEASUREMENTS
96 HOURS TREND PLOT
96 HOURS TREND TABLE

TECHNICAL SPECIFICATIONS

SAFETY INDEX:-

STANDARD ELECTRIC SHOCK RESISTANCE
ANTI-ELECTROSHOCK TYPE CLASS I EQUIPMENT
EMC TYPE CLASS A ANTI-ELECTROSHOCK DEGREE

POWER SPECIFICATIONS:-

AC POWER 100V-240V AC, 50Hz/60Hz ± 1 Hz, 70VA
BATTERY RECHARGEABLE Li-ion BATTERY (11.1 VDC/4000mAh)

APPEARANCE PARAMETERS:-

SIZE (L*H*D) 303mm*160mm*287mm
WEIGHT ABOUT 3.23 Kg

ENVIRONMENTAL CONDITIONS

OPERATION:-

AMBIENT TEMPERATURE	0 C -- +40C
AMBIENT HUMIDITY	20% -- 85%(NO CONDENSATION)
WORKING ALTITUDE	-500 METER -- 4600 METER

SHIPMENT AND STORAGE:-

AMBIENT TEMPERATURE	-20 C -- +60C
AMBIENT HUMIDITY	10% -- 93%
SHIPMENT AND STORAGE ALTITUDE	-500 METER -- 13100 METER



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2, Nadeppally Complex
Kochi-682026
www.meditechcorp.com



PATIENT MONITOR

MTC M319



ABOUT:-

COMPACT AND PORTABLE, ALLOWING-
UNINTERRUPTED MONITORING
CONTINUOUSLY MEASURE SPO₂, NIBP & PR
ANTI-HIGH FREQUENCY ELECTROSURGICAL-
EQUIPMENT

OPERATION:-

AC POWER & BATTERY
LARGE CAPACITY 2000mAh Lithium BATTERY
(6HOURS BACKUP)

DISPLAY:-

4.3" HIGH BRIGHTNESS LCD DISPLAY WITH TRUE-
COLOR WIDE VIEWING CAPACITY

MONITOR:-

BAR GRAPH & PLETHYSMOGRAPH DISPLAY-
INDICATING BLOOD PERFUSION
DUAL MODE DISPLAY CONFIG: BIG FONT MODE &
NORMAL MODE
SUITABLE FOR ADULT, PEDIATRIC & NEONATE
BAR GRAPH WAVEFORM AND FILLED GRAPH-
WAVEFORM

ALARM:-

AUDIO AND VISUAL ALARMS FOR ALL PARAMETERS.
PULSE TONE WITH VARIABLE PITCH TO INDICATE-
RISING OR FALLING OXYGEN SATURATION LEVEL

DATA STORAGE:-

36 HOURS GRAPHICAL AND TABULAR TREND DATA-
STORAGE.

TECHNICAL SPECIFICATIONS

SAFETY INDEX:-

ANTI-ELECTROSHOCK TYPE CLASS I EQUIPMENT
INTERNAL POWERED EQUIPMENT EMC TYPE CLASS A ANTI-
ELECTROSHOCK DEGREE

POWER SPECIFICATIONS:-

AC POWER 100V-240V AC, 50Hz/60Hz $\pm$ 1Hz, 70VA
BATTERY RECHARGEABLE Li-ion BATTERY (11.1 VDC/
2000mAh)

APPEARANCE PARAMETERS:-

SIZE (L*W*H) 208mm*241mm*189mm
WEIGHT 2.4 Kg.

ENVIRONMENTAL CONDITIONS:-

OPERATION

AMBIENT TEMPERATURE 0C -- +40C
AMBIENT HUMIDITY 20% -- 85% (NO CONDENSATION)
ATMOSPHERIC PRESSURE 70kpa -- 106kpa

SHIPMENT AND STORAGE:-

AMBIENT TEMPERATURE -20 C-- +55C
AMBIENT HUMIDITY 10% -- 95%
ATMOSPHERIC PRESSURE 22KPa -- 107.4KPa



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2, Nadeppally Complex
Kochi-682026
www.meditechcorp.com



PATIENT MONITOR

MTC M910



ABOUT:-

SIX PARAMETER (ECG,RESP,SPO2,PR,NIBP,DUAL-TEMP)
BACK LIGHT CONTROL
HIGH FREQUENCY SURGICAL UNIT
DEFIBRILLATION PROTECTION
SUPPORT MASIMO SPO2-AUTHORIZATION-PARTNER
13 KIND OF ARRHYTHMIA ANALYSIS
15 KINDS OF DRUG DOSE CALCULATIONS
TO CONNECT CENTRAL MONITOR & OTHER BED-OBSERVATION

OPERATION:-

AC POWER & BATTERY
LARGE CAPACITY 4 HOURS BATTERY BACKUP

DISPLAY:-

12.1 INCH TFT COLORFUL AND CLEAR DISPLAY

MONITOR:-

DISPLAY MODES: STANDARD, TREND, BIG FONT, OXY CRG, VIEWBED SCREEN, NIBP LIST SCREEN
SUITABLE FOR ADULT, PEDIATRIC, NEONATE (BEDSIDE/ OR/ICU/NICU/PICU)

ALARM:-

THREE LEVEL AUDIBLE AND VISUAL ALARM
DUAL ALARM LIGHT: PHYSIOLOGICAL AND TECHNICAL

DATA STORAGE:-

8000 GROUPS NIBP MEASUREMENTS
680 HOURS TREND DATE & TREND GRAPHS
200 GROUPS ALARM EVENTS REVIEWING
2 HOUR WAVE FORMS REVIEWING

TECHNICAL SPECIFICATIONS

SAFETY INDEX:-

STANDARD ELECTRIC SHOCK RESISTANCE
ANTI-ELECTROSHOCK TYPE CLASS I EQUIPMENT
EMC TYPE CLASS A ANTI- ELECTROSHOCK DEGREE

POWER SPECIFICATIONS:-

AC POWER 100V-240V AC, 50Hz/60Hz, 90VA
BATTERY Li-ion BATTERY (14.8 VDC/2200mAh) 4HOURS
BACKUP

APPEARANCE PARAMETERS:-

SIZE (L*H*D) 300mm*260mm*120mm
WEIGHT ABOUT 3.0 Kg

ENVIRONMENTAL CONDITIONS

OPERATION:-

AMBIENT TEMPERATURE	0 C -- +40C
AMBIENT HUMIDITY	20% -- 85%(NO CONDENSATION)
WORKING ALTITUDE	-500 METER -- 4600
METER	

SHIPMENT AND STORAGE:-

AMBIENT TEMPERATURE	-20 C-- +60C
AMBIENT HUMIDITY	10% -- 93%
SHIPMENT AND STORAGE ALTITUDE	-500 METER -- 13100
METER	



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2, Nadeppally Complex
Kochi-682026
www.meditechcorp.com



SYRINGE PUMP

MTC S213



FEATURES:-

SAFE, CONVINIENT, USER FRIENDLY, EASY TO-
OPERATE
ACCEPT SYRINGES OF 10ml, 20ml, 30ml, AND-
50/60ml
IDENFIES AND DISPLAY THE SYRINGE TYPE-
AUTOMATICALLY
ACCURATE AND SENSITIVE PRESSURE SENSOR-
SHORTENS ALARM RESPONSE TIME AND-
PREVENTS FAST INFUSION DUE TO OCCCLUSION.
VOICE NOTIFICATION FOR ALARM CONDITIONt
DIFFERENT MODES OF OPERATIONS- VOLUME/
TIME, RATE, BODY WEIGHT

OPERATION:-

AC POWER & BATTERY
BATTERY BACKUP : 4 HOURS

DISPLAY:-

LED DISPLAY FOR FLOW RATE
LCD DISPLAY WITH BACKLIGHT FOR MENU-
SETTINGS AND SYSTEM STATUS

ALARMS:-

INFUSION NEAR COMPLETION, INFUSION-
COMPLETE, OCCLUSION, IMPROPER SYRINGE
INSTALLATION, WRONG SETTING, LOW BATTERY,
SYRINGE LOOSE ETC.

TECHNICAL SPECIFICATIONS

APPEARANCE PARAMETERS:-

SIZE (L*W*H) 280mm*210mm*130mm

WEIGHT 2.2 Kg

POWER:-

AC POWER 100V -240 V, 50/60 Hz, 30VA

BATTERY Li-ion BATTERY, 1.6Ah, 4 HOURS BACKUP

OCCCLUSION:-

HIGH 800mmHg $\pm$ 200mmHg (106.7kPa $\pm$ 26.7kPa

MEDIUM 50/600mmHg $\pm$ 100mmHg (66.7kPa $\pm$ 13.3kPa

LOW 300mmHg $\pm$ 100mmHg (40.7kPa $\pm$ 13.3kPa

ALARMS:-

INFUSION NEAR COMPLETION, INFUSION COMPLETE,
OCCLUSION, IMPROPER SYRINGE INSTALLATION, WRONG
SETTING, LOW BATTERY, SYRINGE LOOSE ETC.

RATE:-

50/60ml Syringe : 0.1ml/h – 999.9ml/h (0.1ml/h step)

: 1000ml/h -1500ml/h (1ml/h step)

30ml Syringe : 0.1ml/h – 900ml/h (0.1ml/h step)

20ml Syringe : 0.1ml/h – 600ml/h (0.1ml/h step)

10ml Syringe : 0.1ml/h – 300ml/h (0.1ml/h step)

FLOW RATE ACCURACY WITHIN $\pm$ 3% (AFTER CURRENT CALIBRATION)

MECHANICAL PRECISION WITHIN $\pm$ 2%

TOTAL INJECTION VOLUME 0.1ml – 9999.9ml (0.1ml step)



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2, Nadeppally Complex
Kochi-682026
www.meditechcorp.com



ECG MACHINE

MTC E7100



FEATURES:-

12 LEAD SYNCHRONOUS ACQUISITION AND-DISPLAY
INSTITUTION & PATIENT INFO ENTRY
SAMPLING RATE 32000Hz
BACKLIGHT CONTROL
216mm OR 210mm PRINTOUT
REALTIME ECG FREEZING &REVIEW
PACE DETECTION
HIGH PERFORMANCE CHIP MODULE MAKES-
ACCURATE ECG SIGNAL COLLECTION
CONTINUOUS OPERATION TIME IN BATTERY-
POWER : MORE THAN 8 HOURS

DISPLAY:-

9.7 INCH TILTED HIGH RESOLUTION COLOR-
TOUCH SCREEN LCD DISPLAY

EXTRAS:-

SATDBY MODE &AUTO SHUT DOWN MODE
KEYBOARD, MOUSE, EXTERNAL PRINTER

ECG REPORTS:-

AUTO ECG REPORT WITH &WITHOUT-
INTERPRETATION,
AUTO ECG RHYTHM LEAD ANALYSIS REPORT-
WITH &WITHOUT INTERPRETATION,
MANUAL MODE-REAL TIME RHYTHM.

LEAD LAYOUT OPTION:-

3x4, 3x4+1R, 3x4+3R, 6x2, 6x2+1R, 6x2+3R, 12x1.

SAVING OPTION:-

ECG SAVING OPTION: WITH OR WITHOUT PRINT
INBUILT 1000 ECG STORAGE.
EXPANDABLE MEMORY USB OR SD CARD.
FORMAT: PDF/ PNG/ JPG/ GDT/ DICOM/XML/
HL7/ ZQ ECG

NETWORK:-

DICOM,LAN,WIFI

TECHNICAL SPECIFICATIONS

SAFETY INDEX:-

CLASS I,TYPE CF,WITH PROTECTIVE CIRCUIT OF DEFIBRILLATION-
AND PACE-MAKING

POWER:

AC POWER 100V-240V,50Hz/60Hz ± 1 Hz,120VA
BATTERY RECHARGEABLE Li-ion BATTERY (14.8V/4400mAh)

APPEARANCE PARAMETERS:-

SIZE (L*W*H) 285mm*360mm*142mm
WEIGHT ABOUT 3.9 Kg

ENVIRONMENTAL CONDITIONS

OPERATION:-

AMBIENT TEMPERATURE	+5C -- +40C
AMBIENT HUMIDITY	20% -- 85% (NO CONDENSATION)
ATMOSPHERIC PRESSURE	570hpa – 1060hpa

SHIPMENT AND STORAGE:-

AMBIENT TEMPERATURE	-20 C -- +55C
AMBIENT HUMIDITY	10% -- 95%
ATMOSPHERIC PRESSURE	500hPa – 1060hPa




Contact us now
+91 98467 46764

Meditech Corporations
72/26-B2,Nadeppally Complex
Kochi-682026
www.meditechcorp.com



OUR INSTALLATIONS / CLIENTS IN KERALA



WHY MEDITECH?

- ✓ Competitive Price
- ✓ Turnkey Project Experience
- ✓ Qualified Engineers
- ✓ Largest installations in Kerala
- ✓ Largest installations in Private Hospitals
- ✓ Largest installations in Government Hospitals

MAJOR CLIENTS



CONTACT

MediTech Corporations
72/26-B2, Nadeppally
-Complex, Kochi -682026

www.meditechcorp.com

enquiry@meditechcorp.com

Call Centre: 9846 746 764

