

ECG MACHINE

MTC E617



FEATURES:-

COMPACT & EASY OPERATION
ACCURATE PACE PULSE IDENTIFICATION &
AUTOMATIC ANALYSIS FUNCTION

OPERATION:-

AC POWER & BATTERY
HIGH CAPACITY 4400mAh BATTERY

DISPLAY:-

10 INCH GRAPHIC COLOR DISPLAY (800*480)
12 LEAD SIMULTANEOUS WAVEFORM DISPLAY

ECG REPORTS:-

AUTO/MANUAL/R-R ANALYSIS/STORAGE MODES
AUTOMATIC ANALYSIS FUNCTION
3Ch+++, 6Ch, 6Ch+, 12Ch FORMAT

SAVING OPTION:-

INBUILT 250 ECG STORAGE (REVIEW & PRINT)

NETWORK:-

LAN

TECHNICAL SPECIFICATIONS

SAFETY INDEX:-

CLASS I, TYPE CF, WITH ANTI-DEFIBRILLATION FUNCTION
ELECTROMAGNETIC COMPATIBILITY: CISPR 11, GROUP I, CLASS A

POWER SPECIFICATIONS:-

AC POWER 100V-240V, 50Hz/60Hz $\pm$ 1Hz, 60VA
BATTERY RECHARGEABLE Li-ion BATTERY
(14.4V/2200mAh or 4400mAh)

APPEARANCE PARAMETERS:-

SIZE (L*W*H) 347 mm*293mm*83mm
WEIGHT ABOUT 4.8 Kg

ENVIRONMENTAL CONDITIONS

OPERATION:-

AMBIENT TEMPERATURE	+5C -- +40C
RELATIVE HUMIDITY	25% -- 80% (NO CONDENSATION)
ATMOSPHERIC PRESSURE	700hPa – 1060hPa

SHIPMENT AND STORAGE:-

AMBIENT TEMPERATURE	-20 C-- +55C
AMBIENT HUMIDITY	10% -- 95%
ATMOSPHERIC PRESSURE	500hPa – 1060hPa



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2, Nadeppally Complex
Kochi-682026
www.meditechcorp.com

