

FETAL MONITOR

MTC F372



MAIN FEATURES:-

MAIN BOARD WITH LINUX OS
12-CRYSTAL PULSE DOPPLER TRANSDUCER
EVENT MARKER:-ANNOTATE PERCEIVED FETAL-
MOVEMENTS OR ANY OTHER EVENTS
AUTOMATIC/MANUAL FETAL-
MOVEMENT DETECTION
PATIENT DATA ENTRY OPTION
HEART RATE OFFSET MODE : VISUALLY-
DISTINGUISH BETWEEN TWINS HEART RATE
CROSS CHANNEL VERIFICATION (CCV) FETAL OR-
MATERNAL HEART RATE
NST TIMER

OPERATION:-

AC POWER & BATTERY
LARGE CAPACITY 2200mAh Lithium BATTERY
(4-5 HOURS,WORKING TIME :2 HOURS)

DISPLAY:-

FOLDED 10.4 INCH HIGH BRIGHTNESS TFT LED
DISPLAY

MONITOR:-

FHR1,FHR2,TOCO,AFM WITH WAVEFORM
MODES: STANDARD INTERFACE AND BIG FONT
INTERFACE
FHR SIGNAL QUALITY INDICATOR HELPS OPTIMIZE
THE PROBE POSITION

NETWORK:-

LAN CONNECTIVITY

TECHNICAL SPECIFICATIONS

PRINT:-

FHR1,FHR2,TOCO,AFM
OFFSET OPTION:FHR2 FROM -30 TO +30
PAPER SPEED : 1cm/min, 2cm/min, 3cm/min
CAN SET WAVEFORM WIDTH SIZE:NORMAL,DEMI-BOLD,BOLD
PRINTER : THERMAL 112mm/ Z-fold (REAL TME & REVIEW)

SAVING OPTION:-

8640 HOURS MEMORY FOR FETAL TRACES
PDF EXPORT TO USB

EXTRAS:-

KEYBOARD,MOUSE

SAFETY INDEX:-

EQUIPMENT EMC : IEC60601- CLASS A OF THE GROUP I
CLASS II EQUIPMENT WITH THE ANTI-ELECTRIC SHOCK DEGREE
OF BF

POWER SPECIFICATIONS:-

AC POWER	100V-240V AC,50Hz/60Hz,75VA
ADAPTER	9V DC 5A
BATTERY	RECHARGEABLE Li-ion BATTERY/ 7.4V 2200mAh

APPEARANCE PARAMETERS:-

SIZE	36cm (W) *26cm (D) *8cm (H)
WEIGHT	ABOUT 2.8kg



Contact us now
+91 98467 46764

Meditech Corporations
72/26-B2,Nadeppally Complex
Kochi-682026
www.meditechcorp.com

