

PATIENT MONITOR

MTC M710



ABOUT:-

MEASURED PARAMETER: ECG, RESP, SPO2, PR, NIBP, DUAL TEMP
SIMPLE AND FRIENDLY OPERATING DISPLAY-INTERFACE
PLAYBACK AND BROWSE FUNCTION FOR WAVEFORM AND MONITOR PATIENT DATA
A FULL-SYNCHRONISTIC LEAD MULTICANAL ECG DISPLAY
REAL TIME ST SEGMENT ANALYSIS AND PACEMAKER DETECTION MODE
DYNAMIC WAVEFORMS AND MULTI-LEAD ECG WAVEFORMS DISPLAY
PATIENT INFO INPUT MANAGEMENT FUNCTION
USB AND NETWORK PORT

OPERATION:-

AC POWER & BATTERY
LARGE CAPACITY 4000mAh BATTERY (MORE THAN 60 MINUTES BACKUP)

DISPLAY:-

12.1 INCH WITH TRUE COLOR LCD DISPLAY
WIDE VIEWING ANGLE
HIGH BRIGHTNESS LCD DISPLAY

MONITOR:-

DISPLAY MODES: STANDARD, TREND SCREEN, BIG FONT, OXY CRG
SUITABLE FOR ADULT, PEDIATRIC, NEONATE (BEDSIDE/OR/ICU/NICU/PICU)

DATA STORAGE:-

400 NIBP MEASUREMENTS
96 HOURS TREND PLOT
96 HOURS TREND TABLE

TECHNICAL SPECIFICATIONS

SAFETY INDEX:-

STANDARD ELECTRIC SHOCK RESISTANCE
ANTI-ELECTROSHOCK TYPE CLASS I EQUIPMENT
EMC TYPE CLASS A ANTI-ELECTROSHOCK DEGREE

POWER SPECIFICATIONS:-

AC POWER 100V-240V AC, 50Hz/60Hz ± 1 Hz, 70VA
BATTERY RECHARGEABLE Li-ion BATTERY (11.1 VDC/4000mAh)

APPEARANCE PARAMETERS:-

SIZE (L*H*D) 303mm*160mm*287mm
WEIGHT ABOUT 3.23 Kg

ENVIRONMENTAL CONDITIONS

OPERATION:-

AMBIENT TEMPERATURE	0°C -- +40°C
AMBIENT HUMIDITY	20% -- 85% (NO CONDENSATION)
WORKING ALTITUDE	-500 METER -- 4600 METER

SHIPMENT AND STORAGE:-

AMBIENT TEMPERATURE	-20°C -- +60°C
AMBIENT HUMIDITY	10% -- 93%
SHIPMENT AND STORAGE ALTITUDE	-500 METER -- 13100 METER



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2, Nadeppally Complex
Kochi-682026
www.meditechcorp.com

