

PATIENT MONITOR

MTC M910



ABOUT:-

SIX PARAMETER (ECG,RESP,SPO2,PR,NIBP,DUAL-TEMP)
BACK LIGHT CONTROL
HIGH FREQUENCY SURGICAL UNIT
DEFIBRILLATION PROTECTION
SUPPORT MASIMO SPO2-AUTHORIZATION-PARTNER
13 KIND OF ARRHYTHMIA ANALYSIS
15 KINDS OF DRUG DOSE CALCULATIONS
TO CONNECT CENTRAL MONITOR & OTHER BED-OBSERVATION

OPERATION:-

AC POWER & BATTERY
LARGE CAPACITY 4 HOURS BATTERY BACKUP

DISPLAY:-

12.1 INCH TFT COLORFUL AND CLEAR DISPLAY

MONITOR:-

DISPLAY MODES: STANDARD, TREND, BIG FONT, OXY CRG, VIEWBED SCREEN, NIBP LIST SCREEN
SUITABLE FOR ADULT, PEDIATRIC, NEONATE (BEDSIDE/ OR/ICU/NICU/PICU)

ALARM:-

THREE LEVEL AUDIBLE AND VISUAL ALARM
DUAL ALARM LIGHT: PHYSIOLOGICAL AND TECHNICAL

DATA STORAGE:-

8000 GROUPS NIBP MEASUREMENTS
680 HOURS TREND DATE & TREND GRAPHS
200 GROUPS ALARM EVENTS REVIEWING
2 HOUR WAVE FORMS REVIEWING

TECHNICAL SPECIFICATIONS

SAFETY INDEX:-

STANDARD ELECTRIC SHOCK RESISTANCE
ANTI-ELECTROSHOCK TYPE CLASS I EQUIPMENT
EMC TYPE CLASS A ANTI- ELECTROSHOCK DEGREE

POWER SPECIFICATIONS:-

AC POWER 100V-240V AC, 50Hz/60Hz, 90VA
BATTERY Li-ion BATTERY (14.8 VDC/2200mAh) 4HOURS
BACKUP

APPEARANCE PARAMETERS:-

SIZE (L*H*D) 300mm*260mm*120mm
WEIGHT ABOUT 3.0 Kg

ENVIRONMENTAL CONDITIONS

OPERATION:-

AMBIENT TEMPERATURE	0 C -- +40C
AMBIENT HUMIDITY	20% -- 85%(NO CONDENSATION)
WORKING ALTITUDE	-500 METER -- 4600
METER	

SHIPMENT AND STORAGE:-

AMBIENT TEMPERATURE	-20 C -- +60C
AMBIENT HUMIDITY	10% -- 93%
SHIPMENT AND STORAGE ALTITUDE	-500 METER -- 13100
METER	



Contact us now

+91 98467 46764

Meditech Corporations
72/26-B2, Nadeppally Complex
Kochi-682026
www.meditechcorp.com

